

Beckman Allegra 21R



Beckman Allegra 21R high-speed refrigerated benchtop centrifuge

Rating: Not Rated Yet

Price

Sales price \$3,250.00

[Ask a question about this product](#)

Manufacturer [Beckman](#)

Description

This Beckman Allegra 21R high-speed refrigerated benchtop centrifuge is microprocessor controlled and can get to a top speed of 15,300 rpm w/a fixed angled rotor. Digital display indicates actual speed, set speed, temperature, and rotor imbalance. Accepts many different rotors, which are not included. Very clean.

Please call us at (650) 871-5707 to check availability of this product!