

Beckman Coulter Allegra X-22R Refrigerated Centrifuge



Combining top-notch and uncompromising quality with an affordable price, the Allegra X-22R Centrifuge can fulfill your laboratory needs while still offering a great value for the money.

Rating: Not Rated Yet

Price

Sales price \$4,750.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

Catalog no. 392187, 15,500 rpm maximum speed, 22,065xg maximum g-force, -20°C to +40°C temperature range setable in 1°C increments, features 10 acceleration and deceleration profiles, brushless induction drive, CFC-free refrigeration, variety of rotor options available, 120V, 60Hz, 10.5A.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com