

Beckman Coulter Microfuge 20



Like New! Speed and versatility in a compact footprint.

Rating: Not Rated Yet

Price

Sales price \$1,550.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

Fast, cool and versatile, designed to meet a wide variety of application needs including,

Nucleic acid and protein preparation

Pelleting, extractions

Purifications

Concentrations

Phase separations and receptor binding

Rapid sedimentation of protein precipitates, particulates and cell debris.

15,000 RPM

20,627 x g

120V

60Hz

2.3A