


## Beckman TL-100



Beckman TL-100 benchtop ultracentrifuge to 100,000 rpm. The Classic!!

Rating: Not Rated Yet

**Price**

Sales price \$5,750.00

[Ask a question about this product](#)

Manufacturer [Beckman](#)

### Description

Beckman TL-100 refrigerated ultracentrifuge. Features a microprocessor control, digital display, and vacuum encased induction drive for long reliable performance. Has a max speed of 100,000 RPM and max RCF of 543,000 x g. Temperature range from 2 to 40 degrees C. External dimensions of 29 W x 21 H x 23D

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
(650) 871-5707 office phone & fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)