

## Cole Parmer Masterflex L/S Compact Single-Channel Pump, Model 77240-20 & L/S Miniflex Pump Head



Small footprint; flow rates of 1.6 to 239 mL/min

Rating: Not Rated Yet

### Price

Sales price \$450.00

[Ask a question about this product](#)

Manufacturer [Cole-Parmer](#)

### Description

Quick set-up pump with repeatable low flow rates. Flow rates of 1.6 to 239 mL/min (tubing size selection affects). With power supply and AC cord.

Specs:

Pump Series: L/S Control

Min RPM: 35

Max RPM: 200

Number of Channels: 1 Max Flow Rate Per Channel (330 mL/min):

Tubing sizes accepted: L/S 13, L/S 14, L/S 16, L/S 25

Drive IP rating: IP 22

Number of Rollers: 3

Power (VAC): 90 to 260 Power (Hz): 50/60

Width (in): 5.25; Height (in): 4.5 ; Length (in): 8.5

Width (cm): 13.3 Height (cm): 11.4 Length (cm): 21.6

*Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.*

Biotech Equipment Sales, Inc.

226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotequipmentsales.com](mailto:sales@biotequipmentsales.com)