

Eppendorf 5417R



This, the world's most popular compact refrigerated microcentrifuge, offers variable speed control up to 16,400 rpm (25,000 x g).

Rating: Not Rated Yet

Price

Sales price \$3,600.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

This, the world's most popular compact refrigerated microcentrifuge, offers variable speed control up to 16,400 rpm (25,000 x g). Digital display includes a 99 minute timer. Can also be set for a continuous run. Short spin button for quick spins. A selection of autoclavable rotors for 1.5, 2.0, 0.5, and 0.2 ml tubes is available.

5417 R has a temperature range of -9° to 40°C. A fast cool function cools the centrifuge to 4°C in less than 20 minutes. Measures 12.2" W x 23.5" D x 9.8" H and weighs 75 pounds. 115 V, 60 Hz.