

## Eppendorf 5424 Microcentrifuge



Capacity: 48mL (24 x 1.5/2.0mL) Max. speed/RCF: 15,000rpm/21,130 x g Full digital speed display Variable speed control from 100 to 15,000rpm in 100rpm increments Timer: 30 sec. to 9 hr. 59 min. with HOLD for continuous run Soft br

Rating: Not Rated Yet

### Price

Sales price \$1,550.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

### Description

- Capacity: 48mL (24 x 1.5/2.0mL)
- Max. speed/RCF: 15,000rpm/21,130 x g
- Full digital speed display
- Variable speed control from 100 to 15,000rpm in 100rpm increments
- Timer: 30 sec. to 9 hr. 59 min. with HOLD for continuous run
- Soft brake protects delicate samples
- Short-spin is programmable to desired rpm setting
- Innovative lid design—no sharp edges or exposed parts
- Equipped with 24 x 1.5/2.0mL certified aerosol-tight rotor and lid