

Hettich Kompspin K PC-4 Centrifuge



Table Top K PC-4 Centrifuge adjustable time (up to 99 min) and rotor speed (to 4,500 rpm) with LED displays start, stop and lid release buttons active cycle indicator with power on/off rocker switch on back of centrifuge

Rating: Not Rated Yet

Price

Sales price \$450.00

[Ask a question about this product](#)

Manufacturer [Hettich](#)

Description

Table Top K PC-4 Centrifuge

adjustable time (up to 99 min) and rotor speed (to 4,500 rpm) with LED displays start, stop and lid release buttons active cycle indicator with power on/off rocker switch on back of centrifuge