

## Lancer 815LX Dishwasher



Freestanding dishwasher.

Rating: Not Rated Yet

**Price**

Sales price \$7,000.00

[Ask a question about this product](#)

Manufacturer [Lancer](#)

**Description**

280V, 3 phase. With two racks, as pictured. Undercounter unit, with plinth below to become freestanding.

**Features:**

- Hot, forced air chamber drying
- Efficient direct injection cleaning
- Flexibility with two independent wash levels and interchangeable racks
- Dosing by two peristaltic pumps
- High grade stainless steel clad on all sides for stand alone applications

With plinth that raises the loading level of the 815 LX Glassware Washer by 20 inches (500 mm) and allows for storage of two cleaning chemical

## Miscellaneous: Lancer 815LX Dishwasher

---

containers or rack accessories under the chamber.

### Specifications:

Washer Height: 845 mm / 33.25 in

Washer Width: 600 mm / 23.62 in

Washer Depth: 737 mm / 29 in

Chamber Height: 508 mm / 20 in

Chamber Width: 535 mm / 21 in

Chamber Depth: 510 mm / 20 in

Drying: Hot, forced air drying system

Total electrical power: 7 kw

Net weight: 85 kg / 188 lbs

Loading Height: 250 mm / 10 in

Water Consumption: 12 L / 3.1 gal - per fill

Heat Loss: 2,380 Btu/600 Kcal/h

Sound Level: Steam Condenser: Standard

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
**[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)**