

Vulcon Technologies CS6C Clinaseal Centrifuge



THE IDEAL PHYSICIAN'S OFFICE CENTRIFUGE. The quiet Clinaseal features exclusive AeroLok technology that helps protect personnel from aerosols while maintaining sample temperature stability.

Rating: Not Rated Yet

Price

Sales price \$150.00

[Ask a question about this product](#)

Manufacturer [Vulcon Technologies](#)

Description

THE IDEAL PHYSICIAN'S OFFICE CENTRIFUGE. The quiet Clinaseal features exclusive AeroLok technology that helps protect personnel from aerosols while maintaining sample temperature stability. The motor chamber is isolated from the spinning chamber and the chamber in turn is isolated from the outside environment. A brushless, permanently lubricated motor means no maintenance worries. The sturdy aluminum cabinet is covered with a durable, baked-on powder coating in an attractive tan color that compliments any environment.

Specifications

- Height - 9 in (22.9 cm)
- Width - 9.5 in (24.1 cm)
- Depth - 10 in (25.4 cm)
- Weight - 12 lb (5.4 kg)
- Speed - 3150 RPM
- RCF (relative centrifugal force) - 1132 x g
- Accomodates 2ml-15ml specimen tubes