

AirClean Systems Combination 600 PCR Workstation



Clean-air environment with UV light sterilization. Model # AC632LFUVC

Rating: Not Rated Yet

Price

Sales price \$1,575.00

[Ask a question about this product](#)

Description

SOLD!!!

Features:

UVtect (TM) microprocessor controller

Class 100 (ISO 5) vertical laminar flow air

No joints or gaps in construction

254nm ultraviolet light with user-set digital timer

Constant HEPA-filter monitoring

Interior dimensions (in.): 23D x 31W x 18H