

Coulter Z1



Coulter Z1 Cell and Particle Counter, single threshold, with keyboard and 100 micron aperture tube

Rating: Not Rated Yet

Price

Sales price \$4,850.00

[Ask a question about this product](#)

Manufacturer [Coulter](#)

Description

This Coulter Z1 Cell and Particle Counter, single threshold, with keyboard and 100 micron aperture tube permits the setting of one size threshold. Cells equal to or larger than this setting will be counted. This unit is ideal as a fully automated instrument for the rapid enumeration of both tissue culture and blood cells, as well as obtaining rapid total cell counts on cell lines containing uniformly sized cells. 1-120 micron measurement range. Order starter kit from manufacturer.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com

