

Lindberg/Blue M™ Moldatherm™ Box Furnace 51866A-1

These Box Furnaces include unique insulation and heating element composites.

Rating: Not Rated Yet

Price

Sales price \$1,750.00

[Ask a question about this product](#)

Manufacturer [Lindberg](#)

Description

Minimize outer surface temperatures while maintaining uniform heat distribution within the chamber.

Interior Dimensions (D x W x H): 9 x 6 x 6 in.

Temperature range: 100 - 1100C

Capacity: 5.3 L

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc. is a nine-year-old company based in South San Francisco, the Biotech capital of the world. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.