


New Brunswick Scientific I Series 24R Refrigerating Incubator Shaker



Refrigerating incubator shaker with universal platform.

Rating: Not Rated Yet

Price

Sales price \$4,750.00

[Ask a question about this product](#)

Manufacturer [New Brunswick Scientific](#)

Description

This refrigerating incubating shaker comes with an 18" x 18" universal platform.

Features:

- Microprocessor controls precisely regulate temperature, speed and running time.
- User-friendly keypad with large LED display is clearly visible and easy to set
- Audible and visual alarms alert you to setpoint deviations or end of study
- Wide temperature range: I24: 7 °C above ambient to 60 °C. (I24R: 15 °C below ambient to 60 °C)
- Powerful fan and heater provide rapid heat up and temperature recovery after door opening.
- RS -232 interface for data logging directly to your PC
- Provides the ability to calibrate the speed and temperature directly through the keypad

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.

226 Miller Ave.

South San Francisco, CA 94080

650-871-5707 office

650-276-7487 fax

sales@biotechequipmentsales.com