

Sheldon Laboratory 1408 Economy Vacuum Oven

Durable construction and strongly braced with stainless steel chamber.

Rating: Not Rated Yet

Price

Sales price \$1,350.00

[Ask a question about this product](#)

Manufacturer [Shel Lab](#)

Description

Sheldon Laboratory 1408 Economy Vacuum Oven, 120 Volts, 1.7 cu ft. Capacity

Double plenum design meets CE safety requirements; has a cool outer surface. Cross-flow ventilation; positive latch handles with spring-loaded glass; a small bench top footprint; silicone gaskets.

Temperature range: 15 to 210C

Rise Time: 90 min.

Interior Dimensions: 12 x 20 x 12 in.

Capacity: 1.7 cu ft.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc. is a nine-year-old company based in South San Francisco, the Biotech capital of the world. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.