

A&D GH-202 Analytical Balance



Capacity: 220 g

Rating: Not Rated Yet

Price

Sales price \$2,450.00

[Ask a question about this product](#)

Description

The A&D GH-202 is a high-precision analytical laboratory balance with a glass draft shield designed to protect samples from air currents during weighing. It is suitable for research, analytical testing, and other applications requiring precise mass measurements.

Key specifications

Capacity: 220 g (standard range); 51 g (high-precision range)

Readability: 0.00001 g (0.01 mg) in precision range; 0.0001 g (0.1 mg) in standard range

Calibration: Built-in internal calibration weight; 90 mm diameter weighing pan.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com