

aapptec Sharpfreeze 350 Freeze Dryer



8-place, 15" vacuum chamber; 220V

Rating: Not Rated Yet

Price

Sales price \$15,000.00

[Ask a question about this product](#)

Description

Easy-to-use

No coil

Automatic defrosting

Built-in second trap to prevent damage to the pump

Faster freeze drying
Freeze dry even acidic samples without damage
Process larger volume of samples in a day
On wheels for mobility

Flexible laboratory freeze dryer for peptides, proteins, nucleotides, PNAs and DNAs. It is easy to control through an electronic control panel on the front of the instrument. The 350 model is used for -55 °C temperature. Suitable for lyophilizing HPLC fractions containing acetonitrile or organic solvents.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com