

## ABEC Bioreactor System



**Deprecated:** file\_exists(): Passing null to parameter #1 (\$filename) of type string is deprecated in [/home/beqsaleshost/public\\_html/administrator/components/com\\_virtuemart/helpers/mediahandler.php](/home/beqsaleshost/public_html/administrator/components/com_virtuemart/helpers/mediahandler.php) on line 690



## Fermentors: ABEC Bioreactor System

---

Custom-engineered bioprocessing platform

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

### Description

The ABEC Bioreactor System is a custom-engineered, cGMP-compliant bioprocessing platform from ABEC, a leader in stainless steel, single-use, and hybrid solutions for biopharmaceutical manufacturing. This type of turnkey system is designed for upstream cell culture or microbial fermentation at pilot-to-commercial scale, featuring integrated vessels, piping, agitation, aeration, and advanced automation for precise process control in vaccine, monoclonal antibody, or biologics production. This large-scale stainless steel setup is ideal for high-titer, long-running, or multi-product facilities where durability and performance consistency are paramount.

Inquire for more details; looking for offers. Originally purchased in 2022 for \$14 million. Located in the San Francisco Bay Area.

### Key Features

**High-Capacity Custom Vessels:** Stainless steel bioreactors/fermenters ranging from pilot scale (hundreds of liters) to 30,000 L+ working volumes, optimized with proven geometry, baffles, and high-power impellers for excellent mixing, oxygen transfer (high OTR), and scalability.

**Integrated Automation & Controls:** Modern touchscreen control panels (visible in the photo) with PLC/DCS systems for real-time monitoring and control of critical process parameters (temperature, pH, DO, agitation, gas flow, pressure, etc.), data logging, and recipe management.

**Full Process Integration:** Includes multiple vessels (reactors, media prep/hold tanks), extensive sanitary piping, valves, sensors, and support systems like CIP (Clean-In-Place) skids — enabling seamless upstream-to-downstream workflows in a compact footprint.

**Robust cGMP Design & Compliance:** 316L stainless steel construction, electropolished surfaces, validated for sterility, with features for long-term reliability, easy maintenance, and regulatory compliance (FDA, EMA, etc.). Built for decades of continuous production.

**Turnkey Engineering & Scalability:** Fully customized and factory-tested for the customer's specific process, with ABEC's expertise in seamless scale-up, hybrid integration (if needed), and ongoing support — minimizing validation time and ensuring batch-to-batch reproducibility.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)