

Accucold Pharma-Vac Refrigerator



18.0 cu ft

Rating: Not Rated Yet

Price

Sales price \$1,550.00

[Ask a question about this product](#)

Description

Accucold Pharma-Vac is a purpose-built pharmacy/vaccine refrigerator line meant to support CDC, VFC, and similar storage guidelines.

Model ARG18PV

Key Features / Specifications

Vaccine-range temperature — Adjustable +2 to +8 °C (35.6–46.4 °F). Target stability about ± 1 °C. External microprocessor controller; glycol-buffered product probe plus air sensor; °C/°F display with min/max (15-minute updates); password lock on settings.

Capacity and layout — 18.0 cu ft (510 L) freestanding upright. Cabinet about 78.5 × 27.5 × 26 in (H × W × D); interior ~68.5 × 24.1 × 18.8 in. 7 adjustable plastic-coated wire shelves, interior fan, LED light. Weight ~275 lb. Casters.

Alarms and security — High/low temperature alarms (visual + audible); open-door alarm (typically after ~4 min, adjustable); remote alarm contacts. Factory key lock (2 keys); self-closing reversible door. Optional keypad lock / data logger.

Defrost and refrigerant — Automatic/cycle defrost. Eco hydrocarbon refrigerant (R290 on the 18 cu ft models). 115 V / 60 Hz, ~1.5 A. Hospital-grade cord on many Pharma-Vac units. ETL listed.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com

