

## Across International Jacketed Glass Reactor R50





With Thermo Insulation Jacket, 50 L, 110V, 60 Hz, 140W.

Rating: Not Rated Yet

**Price**

Sales price \$3,750.00

[Ask a question about this product](#)

Description

Across International Jacketed Glass Reactor R50

Power :110V, 60 Hz, 1.7A, 140W

50L capacity

Specifications

Model: R50

Electrical requirements: 110V 60Hz, 140W

Reaction glass vessel capacity: 13 gallons (50L)

Cooling/heating jacket capacity: 4.2 gallons (16L)

Free high quality thermos insulation jacket with open/close viewing window

Drain port ground clearance: 12.5 inches. Unit height (fully assembled): 88 inches

Dimensions: 63" H, 20" L, 35" W

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**

**226 Miller Ave.**

**South San Francisco, CA 94080**

**650-871-5707 office**

**650-276-7487 fax**

[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)