

Alcatel Pascal 2010 SD



SOLD Alcatel Pascal 2010 SD direct-drive vacuum pump

Rating: Not Rated Yet

Price

Sales price \$1,350.00

[Ask a question about this product](#)

Manufacturer [Alcatel](#)

Description

SOLD

This Alcatel Pascal model 2010 SD vacuum pump is direct-drive & 2-stage. It has a pumping speed of 6.8 cfm, and an ultimate total pressure of 2×10^{-3} mbar w/a closed gas ballast. This pump was used for only a single type of application and has been thoroughly tested by our technicians.