


## Amersham Ultrospec 3300 Pro



Amersham Ultrospec 3300 Pro UV/Vis Spectrophotometer

Rating: Not Rated Yet

**Price**

Sales price \$5,850.00

[Ask a question about this product](#)

Manufacturer [Amersham](#)

### Description

**!!!! - RECENTLY SOLD - !!!!**

This Amersham Ultrospec 3300 Pro UV/VIS Spectrophotometer meets the demanding requirements of the Pharmacopoeia in terms of instrument specification and GLP self-diagnostics, with the benefits of being able to be used for rapid, discrete measurements.

It has a bandwidth of  
Just press to read high performance deuterium and tungsten lamps for reduced cost of ownership.

Has built-in applications software with capability for 50 stored methods.

Includes cDNA application software for measuring incorporation of Cy3, Cy5, fluorescein, and user specified dyes into purified microarray hybridization probes and PCR products.

Provides display and print-out information in English, German, French, Italian, Spanish or Russian.