

## AnD GH-202 Analytical Balance



Tall 251mm (9.88") Weighing Chamber

Rating: Not Rated Yet

### Price

Sales price \$2,350.00

[Ask a question about this product](#)

### Description

Provides advanced resolution, is easy to use and calibrate, and offers adjustable settings based on environmental changes.

#### Features:

- Resolution to 0.01mg
- Percentage and Counting Modes, Multiple Weighing Units
- Tall 251mm (9.88") Weighing Chamber
- Built-in Motorized Calibration Weight
- Automatic or Manual Full Digital Calibration
- One Touch Environmental Filters
- GLP/LIMS/ISO Compliant
- Standard Underhook
- Data Memory (up to 200)
- RS-232C and WinCT Software Standard
- Time and Date Technical

## Balances: AnD GH-202 Analytical Balance

---

Specifications:

Capacity: 220g / 51g

Readability: 0.1mg / 0.01mg

Pan Size: 90mm / 3.6" dia.

Percentage Min Div: 0.01%, 0.1% or 1% (auto prompt)

Counting Min Weight: 0.1 mg

Linearity: 0.1 mg / 0.02 mg

Repeatability / Std Dev:  $\pm 0.2$  mg /  $\pm 0.03$  mg

Stabilization Time: Approx. 3.5 sec (0.1 mg) / 8 sec (0.01 mg)

Sensitivity Drift (10°C-30°C):  $\pm 2$ ppm/°C (when automatic self calibration is not used)

Display Refresh: 5 times per second / 10 times per second

Physical Dimensions (W x D x H): 217 x 442 x 316mm / 8.54 x 17.4 x 12.24"

Breeze Break Dimensions (W x D x H): 164 x 148.5 x 251mm / 6.45 x 5.84 x 9.88"

Operating Temperature: 5°C to 40°C (41°F to 104°F), 85% RH or less (no condensation)

Weight: Approx. 8.2 kg / 18 lb

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)