

Baker 6ft. SG604 SterilGARD e3 Biosafety Cabinet



With stand.

Rating: Not Rated Yet

Price

Sales price \$7,000.00

[Ask a question about this product](#)

Manufacturer [Baker](#)

Description

Class II Type A2 (70% air recirculated and 30% exhausted from a common plenum to the room), comes with stand, 115V.

UV light is installed.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com