

## Beckman Airfuge CLS



NEW Beckman Airfuge CLS w/rotor

Rating: Not Rated Yet

**Price**

Sales price \$3,950.00

[Ask a question about this product](#)

Manufacturer [Beckman](#)

### Description

**!!!! - RECENTLY SOLD - !!!!**

This new Beckman Airfuge CLS comes w/a new Chylomicron Rotor and liner. CLS is a compact, air-driven tabletop ultra centrifuge system from Beckman that can get up to 110,000 rpm. It is the first choice for many clinical laboratories and medical professionals who value consistency and accuracy of the results. The CLS is a simple and extremely rapid method of removing chylomicrons (fat particles in lipemic samples) & is perfect for providing accurate results for lipemic clarification. The special Beckman Coulter Chylomicron Rotor and liner can clean up the lipemic samples in 10 minutes, eliminating the need for sendouts. Also can accept any of the other Airfuge rotors

**Please call us at (650) 871-5707 to check availability of this product!**