

Beckman Allegra 6R



Beckman Allegra 6R refrigerated benchtop centrifuge w/GH 3.8 4-place swinging bucket rotor

Rating: Not Rated Yet

Price

Sales price \$4,500.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

Quite possibly the most popular tabletop large capacity centrifuge in the world, and the workhorse of many laboratories, the Beckman Allegra 6R refrigerated benchtop centrifuge is a dependable and adaptable instrument. This Beckman Allegra 6R refrigerated benchtop centrifuge w/GH 3.8 4-place swinging bucket rotor has a capacity of 4 x 750 ml, max speed w/rotor of 3750 rpm, temp range of -5 to 25C, & 30 minute timer.

We stock these centrifuges at all times, along with every rotor, bucket, adapter and microplate carrier you may need for your samples. GH-3.8A autobalancing rotors, GH-3.7 rotors and accessories and all the rest for your Beckman Allegra 6R centrifuge.

[Biotech Equipment Sales, Inc.](#)

[226 Miller Ave.](#)

[South San Francisco, CA 94080](#)

[650-871-5707 office](#)

[650-276-7487 fax](#)

sales@biotechequipmentsales.com