

## Beckman Coulter Allegra 6



Beckman Coulter Allegra 6 High Volume Benchtop Centrifuge with GH-3.8 swinging bucket rotor

Rating: Not Rated Yet

**Price**

Sales price \$3,000.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

### Description

This Beckman Coulter Allegra 6 High Volume Benchtop Centrifuge with GH-3.8 swinging bucket rotor has a capacity of 4x750 ml, can get to 3750 rpm (3200 rcf), has digital rpm readout, brake, and 30 min timer. Clean.

**Please call us at (650) 871-5707 to check availability of this product!**

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
**[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)**