

Beckman Coulter Allegra 6KR



Beckman Coulter Allegra 6KR refrigerating high-capacity kneewell centrifuge w/GH 3.8 swinging bucket rotor & buckets

Rating: Not Rated Yet

Price

Sales price \$4,950.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

This Beckman Coulter Allegra 6KR refrigerating high-capacity kneewell centrifuge w/GH 3.8 swinging bucket rotor & buckets. Variable speed w/digital readout, max speed 3,750 rpm, -10 to 40C temperature range, high, low, and off brake settings. Buckets can accept inserts for many sizes of tubes. Kneewell design takes your benchtop centrifuge off the bench, with easy undercounter storage. Very clean.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com