

## Beckman Coulter Avanti J-E



Beckman Coulter Avanti J-E High Performance refrigerated floor model centrifuge

Rating: Not Rated Yet

**Price**

Sales price \$11,500.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

**!!!SOLD!!!**

This Beckman Coulter Avanti J-E High Performance refrigerated floor model centrifuge is a compact, energy-efficient centrifuge that delivers performance while saving valuable laboratory space. Ideal for cellular, genetic, or protein analysis. Features include a powerful, high-torque switched reluctance (SR) drive for fast separations, two acceleration and three deceleration rates, and a maintenance-free partial vacuum system that reduces friction on the rotor to generate less heat and consume less energy.

Max speed of 21,000 rpm, max capacity of 4 liters, temp range of -10C to 40C, 208/240V, 24 A, 60 Hz. Used, in excellent condition. Rotors not included but many available.

Recently Sold: Beckman Coulter Avanti J-E

---

Please call us at (650) 871-5707 to check availability of this product!

--	--