

Beckman Coulter Avanti JXN-26 Centrifuge, Manufactured 2024





Avanti JXN-26 Centrifuge, 200-240V, 50/60HZ, 24A

Rating: Not Rated Yet

Price

Sales price \$24,750.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

Product No: B34183 High-throughput sample processing. Powerful high speed features, low heat output, imbalance tolerance and quiet performance. Maximum capacity of 6.0L and RPM of 26,000 (82,000 x g). The Friction Reduction System (FRS) for quicker, quieter and cooler runs. Switched Reluctance (SR) drive creates shorter cycle times with high-torque. Stores up-to 1000 user-defined programs.

Rotor sold separately.

Specifications

Platform: Avanti

Maximum Speed: 26,000 rpm

Maximum g-Force: 81770 xg

BioSafe: Yes

Acceleration: 11 acceleration profiles

Deceleration: 12 profiles

Timed Run: Up to 99 hrs, 59 mins

Max Capacity Volume Range: 1 - 6 L

Refrigeration System: Non-CFC, non-ozone depleting refrigerant

Temperature Range(s): From -10°C to +40°C in 1° increments

Rotor ID: Automatic

Zonal Capability: Yes

Pollution Degree: 2

Environmental Conditions: -10°C to 40°C in 1° increments

Installation Category: II

Software: Virtual Network Computing (VNC), Mobile Fuge app for Apple iOS and Android

Material(s): Steel, Plastic

Finish: Powder Coat Paint

Weight: 310 kg(683.4 lb)

Width: 71 cm(28 in)

Depth: 86 cm(33.9 in)

Height: 86 cm(33.9 in)

Clearances: 7.6 cm (sides) & 16 cm (rear)

Humidity Restrictions: 80% (non-condensing)

Power Requirements: 200-240 VAC, 50/60 Hz

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office

650-276-7487 fax
sales@biotechequipmentsales.com