

## Beckman Coulter JA-18 Rotor



Fixed-angle; Accommodates 10 100 mL Tubes

Rating: Not Rated Yet

**Price**

Sales price \$3,250.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

### Description

This ten-place, aluminum, and fixed-angle rotor is for use in Avanti J-26 XP, XPI, J-E, and J6-MI centrifuges. Ideal for high-force, high-volume applications such as pelleting bacteria, cell membranes, and subcellular organelles.

Application: Pelleting Bacteria, Cell Membranes, and Subcellular Organelles

Rotor Accommodates 10 100 mL Tubes

Ten-Place

Material: Aluminum

Fixed-Angle

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

[Biotech Equipment Sales, Inc.](#)

226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)