

## Beckman Coulter JLA-10.500, Never Used



6 x 600 mL; 10,000 rpm; 18,600 x g

Rating: Not Rated Yet

**Price**

Sales price \$7,500.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

### Description

Model: JLA-10.500

Fixed-Angle Rotor, Aluminum

Includes 6 Biosafety Liquid Canisters

Major applications: Large volume for initial processing of bacteria or other cells from fermenters, clearing cell debris from homogenates, or pelleting subcellular organelles and protein precipitates. For use in Avanti J Series and J2 Series centrifuges.

6 x 600 mL

10,000 rpm

18,600 x g

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office

650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)