

Beckman Coulter JLA-8.1000 Fixed-Angle Rotor



Lightly used. 6-place rotor. From 2016.

Rating: Not Rated Yet

Price

Sales price \$7,150.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

Assembly NOTE The JLA-8.1000 rotor can be used in Avanti JXN-26, J-26S XPseries, Avanti J-26 XP series, and in discontinued Avanti J-20 series and Avanti J-20 XP series centrifuges only.

Maximum speed: 8000 RPM

Critical speed range-a: 200 to 400 RPM

Maximum solution density: 1.2 g/mL

Relative Centrifugal Field-b at maximum speed:

...At rmax (222.8 mm): 15,970 g

...At rav (171 mm): 12,300 g

...At rmins (119 mm): 8530 g

k factor at maximum speed (full 1000-ml bottle): 2482

Nominal bottle capacity: 1000 mL

Nominal rotor capacity: 6 liters

Weight of fully loaded rotor (with lid): 34.6 kg (76 lb)

Rotor body material: aluminum

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com