


Beckman Coulter JLA-8.1000 - Never Used!



6 x 1,000 mL; 8,000 rpm; 15, 970 x g

Rating: Not Rated Yet

Price

Sales price \$11,950.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

***** SOLD*****

Model: J-Lite JLA-8.1000 Fixed-Angle Aluminum Rotor

6 x 1,000 mL
8,000 rpm
15, 970 x g

For use in Avanti J-26 XP Series and Avanti J-26S Series centrifuges.

Major applications: General purpose, large-volume processing, pelleting of bacteria, cell organelles, viruses, and precipitates.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com