

Beckman Coulter JS-4.0 4-Place Swinging-Bucket Rotor



4000 rpm; includes the 4 pictured tube adaptors

Rating: Not Rated Yet

Price

Sales price \$4,250.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

This four-place, aluminum swinging bucket rotor is for use with Avanti™ J-26 XP and XPIBK393124 series and J6-MIBK360291 centrifuge.

Ideal for Rapid Sedimentation of Protein Precipitates, Large Particles, Cells, and Cell Debris

Rotor Capacity: Four Liters or Four Blood Bags

Interchangeable with Microplate Carriers

Material: Aluminum

Four-Place Rotors

40000 rpm

Includes the 4 pictured tube adaptors

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com