

## Beckman Coulter Microfuge 22R Refrigerated Centrifuge



SOLD Provides life science researchers with optimized pelleting of DNA, RNA, and protein on the benchtop.

Rating: Not Rated Yet

**Price**

Sales price \$2,500.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

### Description

**SOLD**

A brushless induction drive offers maintenance-free and quiet operation. Standard 24-place fixed angle rotor holds 1.5–2.0mL tubes and features an integral O-ring for safety. Digital display shows set and actual run conditions in time, rpm/rcf, and temperature. Operations in timed, hold, or pulse (short run) can be selected. All interchangeable rotors come standard with aerosol-tight lids, minimizing contamination from broken or leaking tubes. All rotors are autoclavable at 121°C, except for F241.5P ( **BK367187** ), which is autoclavable at 118°C. CFC-free refrigerant.