

Beckman Coulter Type 45 Titanium Fixed Angle Rotor



Type 45 Ti Rotor, Fixed Angle, Titanium, 6x94 mL, 45,000rpm, 235,000xg

Rating: Not Rated Yet

Price

Sales price \$9,500.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

Type 45 Ti Rotor, Fixed Angle, Titanium, 6x94 mL, 45,000rpm, 235,000xg -

Rotor Type: Ultra - Fixed Angle

Maximum Speed: 45,000 rpm

Maximum RCF at r max (g): 235,000 x g

k-Factor: 133

Rotor Material: Titanium

r min: 35.9 mm

r av: 69.8 mm

r max: 103.8 mm

Nominal Rotor Capacity: 564 mL

Number of Rotor Cavities x Nominal Largest Load Capacity: 6 x 94 mL

Nominal Largest Container/Bucket/Carrier/Canister Dimensions: 38 x 102 mm

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com

