


Beckman Coulter Type 70 Ti Rotor



From 2018.

Rating: Not Rated Yet

Price

Sales price \$8,500.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

Fixed-angle titanium rotor.

Max rpm: 70,000

Max g-force: 504,000

No. tubes: 8 x 39ml

Nominal capacity: 312ml

Centrifuge Rotors: Beckman Coulter Type 70 Ti Rotor

r-Min: 39.5

r-Max: 91.9

kfactor: 44

Material: Titanium

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com