

Beckman F0685 Fixed Angle Rotor



Places: 6 Tube Volumes: 85 mL microcentrifuge tubes

Rating: Not Rated Yet

Price

Sales price \$1,100.00

[Ask a question about this product](#)

Manufacturer [Beckman](#)

Description

Major applications: Pelleting subcellular organelles, viruses, bacteria, mitochondria, chloroplasts, or algae.

Places: 6

Tube Volumes: 85 mL microcentrifuge tubes

	Allegra® 64R	Allegra® X-22	Allegra® X-22R
Max RPM	15,500	8,000	10,000
Max RCF, g	26,040	6,953	10,864

Notes: Maximum speed/RCF are different in the Allegra® 64R centrifuge and the Allegra® X-22 centrifuge. This rotor is rated at a lower speed in the Allegra® X-22 and X-22R centrifuges, so *g*-forces and *k* factors are different.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com