

Beckman JA-14 Fixed Angle Rotor with 6 x 250mL



Beckman JA-14 Fixed Angle Rotor with 6 x 250mL

Rating: Not Rated Yet

Price

Sales price \$1,300.00

[Ask a question about this product](#)

Manufacturer [Beckman](#)

Description

Beckman JA-14 Fixed Angle Rotor, Aluminum with Lid, very clean, corrosion-free internals.

Major applications: General-purpose, large-volume and multi-tube processing.

Adapters permit use of up to 30 each 15-mL tubes and up to 72 each 3-mL tubes.

Max. RPM 14,000

Max. g 30,100

k Factor 1,764

Number of Tubes Volume/Size 6 x 250 mL, 62 x 120 mm, 2.5 x 4.8 in

Rotor Capacity 1.5L

For use in Avanti® J Series and J2 Series centrifuges.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com