

Beckman JA14 Fixed Angle Rotor



Fixed angle rotor from Beckman, type JA-14 (SER.E.16782), suitable for Avanti J-serie and J2-serie centrifuges.

Rating: Not Rated Yet

Price

Sales price \$2,250.00

[Ask a question about this product](#)

Manufacturer [Beckman](#)

Description

Details

Used fixed angle rotor from Beckman, type JA-14 (SER.E.16782), suitable for Avanti J-serie and J2-serie centrifuges.
Main applications: General purpose, large volume and multi-tube processing. Adapters permit use of up to 30 tubes.
Material: polished aluminum.
Maximum rotor speed: 14.000 rpm.
Tube dimensions (ø)x(D): 6,2 x 12,0

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com