

## Beckman MLA-130 Fixed-Angle Rotor



MLA-130 Rotor Package, Fixed-Angle Rotor (35 Angle), Titanium. Major applications: Rapid differential sedimentation (pelleting) of small particles such as subcellular organelles and viruses. For use in the Optima MAX-XP Benchtop Ultracentrifuge and Optima MAX Ultracentrifuge.

Rating: Not Rated Yet

### Price

Sales price \$1,500.00

[Ask a question about this product](#)

Manufacturer [Beckman](#)

### Description

Major applications: Rapid differential sedimentation (pelleting) of small particles such as subcellular organelles and viruses. For use in the Optima MAX-XP Benchtop Ultracentrifuge and Optima MAX Ultracentrifuge.

Platform: Tabletop UCF

Rotor Type: Fixed-Angle (Ultracentrifugation)

Maximum RPM: 130,000 rpm

[Biotech Equipment Sales, Inc.](#)  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)