

## Beckman Optima TLX



SOLD Beckman Optima TKX benchtop ultracentrifuge to 120,000 rpm

Rating: Not Rated Yet

**Price**

Sales price \$9,500.00

[Ask a question about this product](#)

Manufacturer [Beckman](#)

### Description

Beckman Optima TLX Ultracentrifuge - **SOLD**

This Beckman benchtop ultracentrifuge is in excellent physical condition.

System features:

- 120,00 rpm capacity
- Adaptable to many rotors
- imbalance-tolerant
- direct vacuum-encased
- air-cooled drives
- rotor imbalance detection
- moisture-purging vacuum pump
- fixed angle rotors with a fluid-containment annulus
- chlorofluorocarbon-free
- solid-state temperature control systems
- automatic rotor identification

- user diagnostics
- programmability with up to five step changes
- ten selectable acceleration/deceleration rates
- microprocessor controlled Instrument Specifications:
  - Temperature control within 2 deg C of set temp.
  - Speed control within 50 rpm of set speed.
- 29w x 23d x 21h in dimensions

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
(650) 871-5707 office phone & fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)