

Biocare HIER Decloaking Chamber



5 discrete temperature settings ranging between 60°C and 110°C with user programmable times

Rating: Not Rated Yet

Price

Sales price \$1,850.00

[Ask a question about this product](#)

Description

Biocare Decloaking Chamber NxGen

The Decloaking Chamber NxGen has been designed for heat-induced epitope retrieval (HIER) and ease of use. It has 5 discrete temperature settings ranging between 60°C and 110°C with user programmable times. The 110°C antigen retrieval protocol can be completed from start to finish in under an hour. With a 72 total slide capacity and only minutes of hands-on time per run, the Decloaking Chamber NxGen offers a walk-away capability similar to fully automated staining instruments. The Decloaking Chamber NxGen transfers run data to a USB drive for export to a user's computer. The run data recorded includes the date and time per run with temperature and pressure readings throughout. With the Decloaking Chamber NxGen recalling the settings from the last run, a quick start of the same protocol is possible.

Designed to optimize and standardize antibody staining procedures and has been engineered to pass strict laboratory safety and quality control requirements. Temperature, pressure and time can be monitored and recorded with the Decloaking Chamber NxGen to produce consistent staining.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com