


BioRad Mini-Protean Tetra System w/ Power Pac HC



The BioRad Mini-Protean Tetra Cell electrophoresis system can run up to 4 mini gels in less than hour or 2-D gels in less than a day.

Rating: Not Rated Yet

Price

Sales price \$625.00

[Ask a question about this product](#)

Manufacturer [BioRad](#)

Description

!!!SOLD!!!

Easy to assemble with a patented sealing mechanism to prevent leaks, the new Mini-PROTEAN Tetra cell system includes a convenient casting stand for hand casting and the ability to run either handcast or precast gels. The instrument's interchangeable module design allows it's tank to be used for protein transfer with a mini trans-blot module.

For added user convenience, the new Tetra cell lid is designed to ensure correct electrode polarity and electrode assembly orientation. Helpful guidelines are added on the tank to indicate how much buffer to add for two-gel and four-gel electrophoresis.

Please call us at (650) 871-5707 to check availability of this product!