

## BioSpec Mini Bead Beater



Shaking Speed: Variable, 2000–3450 strokes/min

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Manufacturer [BioSpec](#)

### Description

The BioSpec Mini-BeadBeater is a high-energy cell disruptor that efficiently homogenizes microbial cells, plant, and animal tissues in sealed microvials using bead-beating technology, ideal for nucleic acid isolation, protein recovery, and dry/cryo-grinding with no cross-contamination.

### Key Specifications:

Capacity: Processes 1–24 screw-cap microvials (0.5, 1.5, 2.0 mL) or one deep-well plate, handling up to 400 mg (wet weight) per vial

Shaking Speed: Variable, 2000–3450 strokes/min, with a 0.75-inch vial displacement for efficient disruption

Timer: Digital, 0–5 minutes with auto-restart, no motor cool-down required between runs

Bead Types: Compatible with glass, zirconia, or steel beads (0.1–2.5 mm) for various sample types, including resistant tissues

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)