

BioTek 405LSRS Microplate Washer



96-well

Rating: Not Rated Yet

Price

Sales price \$6,500.00

[Ask a question about this product](#)

Manufacturer [BioTek](#)

Description

The BioTek 405 LS RS Microplate Washer is the fifth generation BioTek washer and can be used for washing protocols ranging from gentle cell washing to vigorous ELISAs. The 405 LS RS has a simple user interface for creation and execution of wash, dispense, aspirate and maintenance protocols within the extensive onboard software. New functionality, such as built-in automated buffer switching and “quick-change” manifold designs, offer ease of use. The unique patent-pending Ultrasonic Advantage option eliminates the number one cause of assay failure - clogged manifold tubes. Valuable time is saved with this automated maintenance capability.

Model: 405 LSRS

Automates basic dispensing and aspiration

Accommodates 96-well microplates

Low profile, standard height, and deep well plates

Built in ultrasonic cleaner

Specifications:

Precision < 3% CV: 300 μ L/well (96-well washing); < 4% CV: 80 μ L/well (384-well washing)

Plate Type: 96- and 384-well; Low profile and standard height; Solid and filter bottom (filter pore sizes 0.45 to 1.2 μ m) plates

Washing Speed: 96-wells: 300 μ L/well, 3 cycles: < 30 s; 384-wells: 100 μ L/well, 3 cycles: < 80 s; 384-wells: 400 μ L/well, 1 cycle: < 20 s

Type: Microplate Washer

Wash Cycles: 1 to 250

Humidity: 10 to 85%, noncondensing

Microplate Instruments: BioTek 405LSRS Microplate Washer

Includes: Ultrasonic Advantage; 96-well microplate washing
Width (Metric): 35.6 cm
Height (Metric): 25.4 cm
Depth (Metric): 43.2 cm
Voltage: 100/240 V
Bottles: Optional (4 or 10 L)
Dispense Volume: 25 to 3000 μ L/well
Residual Volume: <2 μ L/well (96- & 384-well plates); 96-tube manifold for 96 wells; 192-tube for 384 wells
Soak Time: Programmable in minutes and seconds, up to 60 min.
Weight (Metric): With internal buffer switching - 16.5 kg; Without - 13.5 kg
Flow Rate: High to Low Flow; Optimized rates for cell assays
Frequency: 50/60 Hz
Increments: 1 μ L

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com