


BTX ECM 630 Electro Cell Manipulator w/ Safety Stand



This system incorporates features that make it one of the most technically advanced exponential decay electroporation systems currently available.

Rating: Not Rated Yet

Price

Sales price \$3,750.00

[Ask a question about this product](#)

Manufacturer [BTX](#)

Description

RECENTLY SOLD!!!!

The generator utilizes the BTX Power Platform Technology design and novel digital user interface. The revolutionary Precision Pulse System provides the researcher unparalleled power in controlling the time constant. With the ability to deliver a maximum of 6000 A in the Low Voltage Mode, the ECM 630 is the most powerful generator in its class. Over 200 ECM 600 protocols may be duplicated with this instrument. The additional ECM 630 resistor selection "none" will allow researchers to reproduce protocols from competitive systems lacking resistor settings or reporting "unlimited" resistance. Indeed, the ECM 630 will perform the widest range of electroporation applications among commercial electroporators on the market.

[Biotech Equipment Sales, Inc.](#)

[226 Miller Ave., South San Francisco, CA 94080](#)

[650-871-5707](#)

sales@biotechequipmentsales.com