

## Buchi V-700 Vacuum Pump



Buchi V-700 Vacuum Pump chemically resistant and totally oil free PTFE diaphragm pump.

Rating: Not Rated Yet

**Price**

Sales price \$850.00

[Ask a question about this product](#)

Manufacturer [Buchi](#)

### Description

**!!!! - RECENTLY SOLD - !!!!**

This Buchi V-700 Vacuum Pump chemically resistant and totally oil free PTFE diaphragm pump designed for the use with all gases and vapours typically associated with the laboratory. With a flow rate of 1.8 cubic meter/h and an ultimate vacuum of below 10 mbar, it covers a wide range of classical laboratory applications and is optimally suited for work with a rotary evaporator. Other typical applications include the evacuation of drying ovens or vacuum centrifuges and the generation of vacuum for filtration purposes.

**Please call us at (650) 871-5707 to check availability of this product!**