

## Carbolite ELF 11/23 Furnace to 1100C



Heat-Up Time: Approximately 49 minutes to 1000°C

Rating: Not Rated Yet

**Price**

Sales price \$2,650.00

[Ask a question about this product](#)

### Description

The Carbolite ELF 11/23 is a bench-mounted laboratory furnace designed for light-duty heat treatment up to 1100°C, featuring a 23-liter chamber, low thermal mass insulation, and a drop-down door for efficient heating and easy sample loading.

### Key Specifications

Electrical requirements: 208V, 50/60 Hz, 5000W, 25A, 23L

Overall Depth, Height, Width: 27 1/4 in, 28 in, 19 7/8 in

Convection Type: Gravity

Capacity: 0.81 cu ft; 23 L

Maximum Temperature: 1100°C

## Ovens/Furnaces: Carbolite ELF 11/23 Furnace to 1100C

---

Chamber Volume: 23 liters

Heat-Up Time: Approximately 49 minutes to 1000°C

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)